

# CHARLOTTE™

SHOULDER SYSTEM




**ARTHROTEK®**  
A BIOMET COMPANY



# Finally...

# a complete


# INSTRUMENTS



**Tray 1**  
Hand Instruments  
and RC Retractor




**Tray 2**  
Portal Placement/  
Miscellaneous Tray



**Tray 3**  
Anchor Instrumentation/  
Suture Management

# shoulder system.



## IMPLANTS

### CHARLOTTE™ SHOULDER SYSTEM

In the last decade, orthopaedic surgeons have witnessed the development of arthroscopic techniques to manage labral tears, shoulder instability, and rotator cuff tears. Clinical results are now being reported that compare favorably to our traditional open procedures. As we learn and adapt these arthroscopic techniques to our own practices, operating rooms must have the necessary instrumentation available to efficiently address these shoulder lesions.

We have all experienced sending the circulating nurse scurrying out of the room to search for the equipment needed to repair an unanticipated labral lesion or to accomplish an arthroscopic rotator cuff repair. The Charlotte™ Shoulder System was designed to eliminate this problem via a comprehensive and versatile shoulder instrument set which can be utilized for carrying out the spectrum of arthroscopic and mini-open shoulder procedures. Resorbable screw-in suture anchors (for both cortical and cancellous bone) and disposable cannulas are available to complement the system.

We hope this system assists you and your operating room staff in performing successful arthroscopic shoulder surgery.

**Patrick M. Connor, M.D.**

**Donald F. D'Alessandro, M.D.**

*Ortho Carolina, Charlotte, NC*



### Precision Hand Instruments

- Over 50 different instruments to choose from
- Tapered tube and shaft
- Oversized, ergonomic handle
- Dual cutting tips

#### **Biters**

- Low profile, biting basket forceps allow easy access in tight spaces
- Many different types and angles to choose from including BackBiter<sup>®</sup> device right and left biter

#### **ArthroPasser™ Suture Retriever**

- Three different angles: 0°, 15° up, 45° up
- passes any type of suture
- sharp tip to assist in passing suture through tissue

#### **Suture Grasper**

- Dual-purpose ring suture grasper can either grasp the suture directly or retrieve suture by allowing it to slide freely through the ring

# CHARLOTTE™

## SHOULDER SYSTEM



#### **Caspari™ Suture Punch**


- Ability to pass any type of suture through soft tissue enabling placement of either mattress or simple stitches
- Does not require the use of shuttle relays
- Two sizes and various angles
- Can be utilized with any anchor that has a sliding suture

#### **Nordt™ Knot Tightener**

- Allows both pushing and spreading of knots, reproducing the same technique used when tying traditional open knots.
- It can be used with any type of suture, and obviates the need to thread the suture through an eyelet








### RC Retractor


- Several blade sizes allow for superficial or deep muscle retraction.
- Permits retraction in three different directions (medial, lateral, and inferior), ideal for open and mini-open repairs.
- The fiberoptic light cable is compatible with all arthroscopic light sources.



### MaxCutter™ Suture Cutter


- Cuts any size and type of suture
- Ergonomic spring loaded handle
- Built-in knot pusher





### Guide Pin

- The guide pin is used to establish the location and direction of unique portals (eg, lateral portal utilized for SLAP lesion repair, anterior trans-subscapularis portal, accessory posterior portal, etc). The cannulated obturators are then placed over the guide pin for controlled cannula placement
- The guide pin has both sharp and blunt ends, with an easily utilized locking handle.





### Calibrated Probe

- Ergonomically designed calibrated probe allows the surgeon to accurately quantify pathology.

### Acromial Rasp

- Acromial rasp has aggressive teeth on one side and finer teeth on the other. In addition, offset and straight ends facilitate smooth, parallel planing of bony resection.




### **LactoScrew® Instruments**

- Drill-tap for all LactoScrew® Anchor sizes
- Ratcheting handle
- Obturator and guide for labral tears
- 3.5 mm drill for extremely dense bone
- 5.5mm awl for dense bone (to be used in conjunction with 5.5mm tap)



# LactoScrew® Suture Anchors

- 
- 2.8, 3.5, 5.5, 6.8mm sizes for various indications
  - LactoSorb® L15 material
  - Simple insertion—Drill-tap only
  - 5.5, 6.8mm available with or without needles
  - Available with MaxBraid™ PE suture, Osteoprene™ resorbable suture or polyester suture

## LactoSorb® L15 Material

- Low crystalline material
- Amorphous
- Maintains 80% of its strength at eight months<sup>1</sup>
- Resorbs in approximately 16–24 months<sup>1</sup>
- Clear radiographs
- Does not retard pediatric skeletal growth

1. Data on file at Arthrotek, Inc.



# Ti-Screw Suture Anchors

## Ti-Screw Suture Anchors with EasySlide™ Process

- 3.0, 5.0, 6.5mm sizes for various indications
- One step insertion – No drill or tap
- Can be repositioned if unsatisfied with initial location
- Tremendous pull-out strength<sup>1</sup>
- All sizes are available with two sutures
- All sizes are available with or without needles
- Available with MaxBraid™ PE Suture or polyester Suture

## Ti-Screw SP Suture Anchors with EasySlide™ Process

- Independent eyelets ease suture sliding and management
- Double or triple suture loaded
- Fully threaded to achieve cortical fixation
- 5.0mm and 6.5mm options

## EasySlide™ Surface Treatment

- More than doubles (132% increase) suture life during sliding<sup>1</sup>
- Increased performance of sliding knots
- Minimizes suture abrasion
- Does not alter the chemical or biological characteristics of the titanium alloy.

1. Data on file at Arthrotek, Inc.

2. Data supplied by DSM Dyneema

3. Data on file at Arthrotek, Inc.



## MaxBraid™ PE Suture

- 100% polyethylene
- Braided Dyneema® Purity Fibers are 15x stronger than steel on a weight for weight basis<sup>2</sup>
- Silky smooth feel is gentle on tissue and gloves
- Unique braid allows suture to lie flat and improves knot security



## Osteoprene™ Resorbable Suture

- One of two sutures indicated for orthopedic indications
- Composed of 88% L-Lactide and 12% Trimethylene Carbonate (TMC)
- 73% loss of strength at three months and 52% at six months<sup>3</sup>
- Available as a #2 in white or violet
- Available on 3.5, 5.5mm LactoScrew® Suture Anchors



# Package Inserts

## Arthrotek, Inc.

A Wholly Owned Subsidiary of Biomet, Inc.  
56 East Bell Drive  
Warsaw, Indiana 46582 USA

01-50-1072

Date: 02/05

### ARTHROTEK® SOFT TISSUE ANCHORING DEVICES ATTENTION OPERATING SURGEON

#### DESCRIPTION

The Arthrotek® Soft Tissue Anchoring Devices are resorbable repair devices used to attach soft tissue to bone. LactoSorb® Soft Tissue Screw and Washer and MicroMax™ Suture Anchor are used with or without a suture. LactoSorb® L-15 Screw Anchors consist of a screw and head design and are used with a suture. The devices are implanted into a predrilled bone hole and are made of a resorbable copolymer, a polyester derivative of lactic acid and glycolic acid. Poly(lactic/polyglycolic acid) copolymer degrades and resorbs *in vivo* by hydrolysis to lactic and glycolic acids, which are then metabolized by the body.

#### MATERIALS

Poly-L-Lactic Acid/Polyglycolic Acid  
Ultra-High Molecular Weight Polyethylene (UHMWPE)  
Polyester  
Polypropylene

#### INDICATIONS

1. LactoScrew® L-15 Screw Anchor (85% PLLA/15% PGA);

##### Shoulder

Bankart repair  
SLAP lesion repair  
Acromio-clavicular separation  
Rotator cuff repair  
Capsule repair or capsulolabral reconstruction  
Biceps tenodesis  
Deltoid repair

##### Wrist/Hand

Scapholunate ligament reconstruction  
Ulnar/radial collateral ligament reconstruction

##### Ankle/Foot

Lateral stabilization  
Medial stabilization  
Achilles tendon repair/reconstruction  
Hallux valgus reconstruction  
Mid- and forefoot reconstruction

##### Elbow

Tennis elbow repair  
Ulnar or radial collateral ligament reconstruction  
Biceps tendon reconstruction

##### Knee

Medial collateral ligament repair  
Lateral collateral ligament repair  
Posterior oblique repair  
Joint capsule closure  
Iliotibial band tenodesis  
Patellar ligament/tendon repair

2. LactoSorb® L-15 Screw and Washer (85% PLLA/15% PGA) and MicroMax™ Suture Anchor:

##### Shoulder Indications

Bankart Repair  
SLAP Lesion Repair  
Acromio-clavicular Separation Repair  
Rotator Cuff Repair  
Capsule Repair and Capsulolabral Reconstruction  
Biceps Tenodesis  
Deltoid Repair

##### Wrist Indications

Scapholunate ligament reconstruction

##### Elbow Indications

Tennis Elbow Repair  
Biceps Tendon Reattachment  
Medial and Lateral Repairs  
Ulnar or Radial Collateral Ligament Reconstruction

##### Knee Indications

Extra-Capsular Repair  
Medial Collateral Ligament Repair  
Lateral Collateral Ligament Repair  
Posterior Oblique Ligament Repair  
Joint Capsule Closure  
Iliotibial Band Tenodesis Reconstruction  
Patellar Ligament/Tendon Repair  
Vastus Medialis Obliquus (VMO) Muscle Advancement

LactoSorb® L-15 Screw and Washer and MicroMax™ Suture Anchor are preloaded with suture for use at the discretion of the physician.

#### CONTRAINDICATIONS

1. Active infection.
2. Patients with mental or neurologic conditions who are unwilling or incapable of following postoperative care instructions.
3. Patient conditions including: blood supply limitations, insufficient quantity or quality of bone, or latent infections.
4. Pathologic soft tissue conditions, which would prevent secure fixation.

#### WARNINGS

Arthrotek® internal fixation devices provide the surgeon with a means to aid in the management of soft tissue to bone reattachment procedures. While these devices are generally successful in attaining these goals, they cannot be expected to replace normal healthy soft tissue or withstand the stress placed upon the device by full or partial weight bearing or load bearing, particularly in the presence of incomplete healing. Therefore, it is important that immobilization (use of external support, sling, etc.) of the treatment site be maintained until healing has occurred. Surgical implants are subject to repeated stresses in use, which can result in fracture or damage to the implant. Factors such as the patient's activity level and adherence to weight bearing or load bearing instructions have an effect on the service life of the implant. The surgeon must be thoroughly knowledgeable not only in the medical and surgical aspects of the implant, but also must be aware of the mechanical and polymeric aspects of the surgical implants.

1. Correct selection of the implant is extremely important. The potential for success in soft tissue to bone fixation is increased by the selection of the proper type of implant. While proper selection can help minimize risks, the device is not designed to withstand the unsupported stress of full weight bearing, load bearing or excessive activity.
2. The implants can loosen or be damaged when subjected to increased loading associated with inadequate healing. If healing is delayed, or does not occur, the implant or the procedure may fail. Loads produced by weight bearing and activity levels may dictate the longevity of the implant.
3. Inadequate fixation at the time of surgery can increase the risk of loosening and migration of the device or tissue supported by the device. Sufficient bone quantity and quality are important for adequate fixation and success of the procedure. Bone quality must be assessed at the time of surgery. Adequate fixation in diseased bone may be more difficult. Patients with poor quality bone, such as osteoporotic bone, are at greater risk of device loosening and procedure failure.
4. Care is to be taken to assure adequate soft tissue fixation at the time of surgery. Failure to achieve adequate fixation or improper positioning or placement of the device can contribute to a subsequent undesirable result.
5. The use of appropriate immobilization and postoperative management is indicated as part of the treatment until healing has occurred.
6. Correct handling of implants is extremely important. Do not modify implants. Do not notch or bend implants. Notches or scratches put in the implant during the course of surgery may contribute to breakage. Intraoperative fracture of devices can occur if excessive force (torque) is applied while seating.
7. DO NOT USE if there is loss of sterility of the device.
8. Discard and DO NOT USE opened or damaged devices, and use only devices that are packaged in unopened or undamaged containers.
9. Ensure contact of tissue to bone when implanting. DO NOT OVERTIGHTEN the screw. Structural damage to the tissue and implant may occur if the screw is overtightened.
10. Adequately instruct the patient. Postoperative care is important. The patient's ability and willingness to follow instructions is one of the most important aspects of successful soft tissue management. Patients affected with senility, mental illness, alcoholism, and drug abuse may be at a higher risk of device or procedure failure. These patients may ignore instructions and activity restrictions. The patient is to be instructed in the use of external supports that are intended to immobilize the repair site

and limit weight bearing or load bearing. The patient is to be made fully aware and warned that the device does not replace normal healthy tissue, and that the device can break, bend or be damaged as a result of stress, activity, load bearing, or weight bearing. The patient is to be made aware and warned of general surgical risks, possible adverse effects, and to follow the instructions of the treating physician. The patient is to be advised of the need for regular postoperative follow-up examination as long as the device remains implanted.

11. MicroMax™ Suture Anchor – Loss of bone fixation may occur if flanged wings are not properly deployed.

#### PRECAUTIONS

Instruments are available to aid in the accurate implantation of internal fixation devices. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments, which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Arthrotek recommends that all instruments be regularly inspected for wear and disfigurement.

If device contains MaxBraid™ suture, refer to manufacturer package insert for further information.

#### POSSIBLE ADVERSE EFFECTS

1. Infection can lead to failure of the procedure.
2. Neurovascular injuries can occur due to surgical trauma.
3. Bending, fracture, loosening, rubbing, and migration of the implant may occur as a result of excessive activity, trauma, or load bearing.
4. Implantation of foreign materials can result in an inflammatory response or allergic reaction.
5. Inadequate healing.
6. Pain, discomfort, or abnormal sensation due to the presence of the device.
7. Necrosis of the bone or tissue.

#### STERILITY

Arthrotek® resorbable implants are sterilized by exposure to Ethylene Oxide (ETO) Gas. Do not resterilize. Do not use past expiration date.

**STORE AT OR BELOW ROOM TEMPERATURE. DO NOT EXPOSE PRODUCT TO TEMPERATURES GREATER THAN 120°F OR 49°C.**

**Caution:** Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

Authorized Representative:  
Biomet U.K., Ltd.  
Waterton Industrial Estates  
Bridgend, South Wales  
CF31 3XA U.K.  
☎ 0086

The information contained in this package insert was current on the date this brochure was printed. However, the package insert may have been revised after that date. To obtain a current package insert, please contact Arthrotek at the contact information provided herein.

**Distributed by:****Arthrotek, Inc.**

A Wholly Owned Subsidiary of Biomet, Inc.

56 East Bell Drive

Warsaw, Indiana 46582 USA

**01-50-1134****Date: 08/05****MaxBraid™**

POLYETHYLENE SUTURE

Non-absorbable Surgical Suture  
U.S.P. except for oversized diameter.**Sterile:** Contents sterile unless package has been opened or damaged. Single Use Only, Do Not Resterilize.**DESCRIPTION**

MaxBraid™ Polyethylene Surgical Suture is a nonabsorbable, sterile surgical suture composed of ultra high molecular weight polyethylene. MaxBraid™ Polyethylene Suture is provided braided as a co-braid of undyed polyethylene and blue monofilament polypropylene or undyed (white). MaxBraid™ sutures are U.S.P. except for diameter in the following sizes:

Suture Size	USP Ave Diameter Specification (mm) <861>	Maximum Oversize Average Diameter (mm)	Maximum Oversize Average Diameter from USP (mm)
2-0	0.300 – 0.339	0.363	0.024
0	0.350 – 0.399	0.459	0.060
1	0.400 – 0.499	0.574	0.075
2	0.500 – 0.599	0.617	0.018

**MaxBraid™ sutures exceed USP specifications for diameter.****INDICATIONS**

MaxBraid™ Polyethylene Nonabsorbable Surgical Sutures are indicated for use in approximation and/or ligation of soft tissues, including use of allograft tissue for orthopedic surgeries.

**ACTIONS**

MaxBraid™ Polyethylene Surgical Suture elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

**WARNINGS**

**Do not resterilize. Do not use if package is opened or damaged. Discard open, unused sutures.**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Polyethylene Surgical Suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

**PRECAUTIONS**

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

**ADVERSE REACTIONS**

Adverse effects associated with the use of this device include: wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

**HOW SUPPLIED**

MaxBraid™ Polyethylene Surgical Suture is available in USP sizes 2-0 through 2 (metric sizes 3 through 5). MaxBraid™ Surgical Suture is provided sterile as a co-braid of undyed polyethylene and blue monofilament polypropylene or undyed (white).

The suture is provided in a variety of lengths, with and without needles, with or without pledgets, and may be supplied in a variety of cut lengths or on ligating reels. Finished suture may be packaged in cartons as single packs, multipacks, or procedure packs.

**STERILITY**

MaxBraid™ Polyethylene Nonabsorbable Surgical Sutures are sterilized by exposure to Ethylene Oxide (ETO) Gas. Do not resterilize. Do not use past expiration date.

MaxBraid™ is a trademark of Arthrotek, Inc.

**Caution:** Federal law (USA) limits this device to sale or use by or on the order of a physician.

Comments regarding this device can be directed to Attn: Regulatory Dept., Biomet Inc., P.O. Box 587, Warsaw, IN 46581 USA, Fax: 574-372-1683.

Manufacturer: Teleflex Medical  
600 Airport Road  
Fall River, MA 02720 USA  
508-677-6600

Telephone 800-367-7874 (USA only)  
+1-508-677-6600

Suture CE marked by Teleflex

**The information contained in this package insert was current on the date this brochure was printed. However, the package insert may have been revised after that date. To obtain a current package insert, please contact Arthrotek at the contact information provided herein.**

# Package Inserts

Arthrotek, Inc. 01-50-1018  
A Wholly Owned Subsidiary of Biomet, Inc.  
P.O. Box 587 Date: 09/05  
56 East Bell Drive  
Warsaw, Indiana 46581 USA

## Arthrotek® Internal Fixation Devices Attention Operating Surgeon

### DESCRIPTION

Arthrotek manufactures a variety of internal fixation devices intended to aid in arthroscopic and orthopedic reconstructive procedures requiring soft tissue fixation, due to injury or degenerative disease. Implants used for this application include: screws, washers, anchors, pins, and suture. Specialty implants are available for specialized treatments.

### Materials:

316 LVM Stainless Steel  
Titanium Alloy  
Ultra-High Molecular Weight Polyethylene (UHMWPE)  
Polyester

### INDICATIONS

The Metal Screw Anchor and the Harpoon® Suture Anchor are indicated for use in soft tissue reattachment procedures in the shoulder, wrist, elbow, and knee. Specific indications as follows:

**Shoulder Indications**—Bankart repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule repair or capsulolabral reconstruction, biceps tenodesis, deltoid repair.

**Wrist/Hand Indications**—Scapholunate ligament reconstruction, ulnar/radial collateral ligament reconstruction.

**Ankle/Foot Indications**—Lateral stabilization, medial stabilization, Achilles tendon repair/reconstruction, hallux valgus reconstruction, mid- and forefoot reconstruction.

**Elbow Indications**—Ulnar or radial collateral ligament reconstruction, biceps tendon reconstruction.

**Knee Indications**—Medial collateral ligament repair, lateral collateral ligament repair, posterior oblique ligament repair, joint capsule closure, iliotibial band tenodesis, and patellar ligament/tendon repair.

Bone Mulch® Screws are intended for use in fixation of semitendinosus and/or gracile tendon grafts in ACL reconstruction, only.

Interference Screws and Set Screws are intended for use in fixation of patellar bone-tendon-bone grafts in ACL reconstruction.

Screw and Washers are indicated for soft tissue fixation to bone, and bone to bone fixation in orthopedic procedures specifically during Ligament reconstruction.

Toggle anchors (ie. toggle buttons and EZLoc™) are indicated for use for fixation of tendons and ligaments during orthopedic reconstruction procedures such as Anterior Cruciate Ligament (ACL) Reconstruction.

Patient selection factors to be considered include: 1) need for soft tissue to bone fixation, 2) ability and willingness of the patient to follow postoperative care instructions until healing is complete, and 3) a good nutritional state of the patient.

### CONTRAINDICATIONS

1. Infection.
2. Patient conditions including blood supply limitations, and insufficient quantity or quality of bone or soft tissue.
3. Patients with mental or neurologic conditions who are unwilling or incapable of following postoperative care instructions.
4. Foreign body sensitivity. Where material sensitivity is suspected, testing is to be completed prior to implantation of the device.

### WARNINGS

Arthrotek® internal fixation devices provide the surgeon with a means to aid in the management of soft tissue to bone reattachment procedures. While these devices are generally successful in attaining these goals, they cannot be expected to replace normal healthy bone or withstand the stress placed upon the device by full or partial weight bearing or load bearing, particularly in the presence of nonunion, delayed union, or incomplete healing. Therefore, it is important that immobilization (use of external support, walking aids, braces, etc.) of the treatment site be maintained until healing has occurred. Surgical implants are subject to repeated stresses in use, which can result in fracture or damage to the implant. Factors such as the patient's weight, activity level, and adherence to weight bearing or load bearing instructions have an effect on the service life of the implant. The surgeon must be thoroughly knowledgeable not only in the medical and surgical aspects of the implant, but also must be aware of the mechanical and metallurgical aspects of the surgical implants.

1. Correct selection of the implant is extremely important. The potential for success in soft tissue to bone fixation is increased by the selection of the proper type of implant. While proper selection can help minimize risks, neither the device nor grafts, when used are designed to withstand the unsupported stress of full weight bearing, load bearing or excessive activity.
2. The implants can loosen or be damaged and the graft can fail when subjected to increased loading associated with nonunion or delayed union. If healing is delayed, or does not occur, the implant or the procedure may fail. Loads produced by weight bearing, and activity levels may dictate the longevity of the implant.
3. Inadequate fixation at the time of surgery can increase the risk of loosening and migration of the device or tissue supported by the device. Sufficient bone quantity and quality are important to adequate fixation and success of the procedure. Bone quality must be assessed at the time of surgery. Adequate fixation in diseased bone may be more difficult. Patients with poor quality bone, such as osteoporotic bone, are at greater risk of device loosening and procedure failure.
4. Implant materials are subject to corrosion. Implanting metals and alloys subjects them to constant changing environments of salts, acids, and alkalis that can cause corrosion. Putting dissimilar metals and alloys in contact with each other can accelerate the corrosion process that may enhance fracture of implants. Every effort should be made to use compatible metals and alloys when marrying them to a common goal, i.e., screws and plates.
5. Correct handling of implants is extremely important. Do not modify implants. Do not notch or bend implants. Notches or scratches put in the implant during the course of surgery may contribute to breakage. Intraoperative fracture of screws can occur if excessive force (torque) is applied while seating bone screws.
6. Do not use excessive force when inserting suture anchors. Excessive force (long hard hammer blows) may cause fracture or bending of the device. When encountering hard cortical bone, predrill with a 3/32 or 1/8 inch drill prior to inserting suture anchors.
7. Adequately instruct the patient. Postoperative care is important. The patient's ability and willingness to follow instructions is one of the most important aspects of successful fracture management. Patients effected with senility, mental illness, alcoholism, and drug abuse may be at a higher risk of device or procedure failure. These patients may ignore instructions and activity restrictions. The patient is to be instructed in the use of external supports, walking aids, and braces that are intended to immobilize the fracture site and limit weight bearing or load bearing.

The patient is to be made fully aware and warned that the device does not replace normal healthy bone, and that the device can break, bend or be damaged as a result of stress, activity, load bearing, or weight bearing. The patient is to be made aware and warned of general surgical risks, possible adverse effects, and to follow the instructions of the treating physician. The patient is to be advised of the need for regular postoperative follow-up examination as long as the device remains implanted.

### PRECAUTIONS

Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat with implants that have been even momentarily placed in a different patient.

Instruments are available to aid in the accurate implantation of internal fixation devices. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments, which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Arthrotek recommends that all instruments be regularly inspected for wear and disfigurement.

If device contains MaxBraid™ suture, refer to manufacturer package insert for further information.

### POSSIBLE ADVERSE EFFECTS

1. Nonunion or delayed union, which may lead to breakage of the implant.
2. Bending or fracture of the implant.
3. Loosening or migration of the implant.
4. Metal sensitivity, or allergic reaction to a foreign body.
5. Pain, discomfort, or abnormal sensation due to the presence of the device.
6. Nerve damage due to surgical trauma.
7. Necrosis of bone or tissue.
8. Inadequate healing.
9. Intraoperative or postoperative bone fracture and/or postoperative pain.

### STERILITY

Arthrotek® internal fixation implants are typically supplied sterile and are sterilized by exposure to a minimum dose of 25kGy of gamma radiation or by Ethylene Oxide Gas (ETO) if device contains MaxBraid™ PE suture. If supplied sterile, do not resterilize the implant.

If not supplied sterile, metallic internal fixation devices must be sterilized prior to surgical use. Do not sterilize UHMWPE implants using steam autoclaving methods. Do not use implants after expiration date.

Pre-VacuumSteam (HI-VAC) -- wrapped or unwrapped  
Temperature 270°-275° F (132°-135°C)  
Exposure Time 5 Minutes  
Drying Time 8 Minutes

Since Arthrotek is not familiar with individual hospital handling methods, cleaning methods and bioburden, Arthrotek cannot assume responsibility for sterility even though the guideline is followed.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Authorized Representative: Biomet U.K., Ltd.  
Waterton Industrial Estates,  
Bridgend, South Wales  
CF31 3XA, U.K.

CE0086

The information contained in this package insert was current on the date this brochure was printed. However, the package insert may have been revised after that date. To obtain a current package insert, please contact Arthrotek at the contact information provided herein.



# Ordering Information

## 2.8mm LactoScrew® Suture Anchors

<b>905494</b>	Long with one #2 MaxBraid™ Suture
<b>905495</b>	Short with One #2 MaxBraid™ Suture with Needles

## 3.5mm LactoScrew® Suture Anchors

<b>905573</b>	One #2 Blue Polyester Suture
<b>905570</b>	One #2 MaxBraid™ Suture
<b>905591</b>	One #2 MaxBraid™ Suture with AutoKnot™ Pre-Tied Knot

## 3.5mm LactoScrew® LS Suture Anchors

<b>905482</b>	Two 2-0 MaxBraid™ Sutures
<b>905493</b>	One, #2 MaxBraid™ Suture
<b>905483</b>	One #2 Polyester Suture
<b>CP9055593</b>	One #2 Osteoprene™ Resorbable Suture

## 5.5mm LactoScrew® Suture Anchors

<b>905571</b>	Two #2 MaxBraid™ Sutures with Needles
<b>905572</b>	Two #2 MaxBraid™ Sutures
<b>905575</b>	Two #2 MaxBraid™ Sutures with Cutting Needles
<b>905592</b>	Two #2 MaxBraid™ Sutures with AutoKnot™ Pre-Tied Knot
<b>905491</b>	Two #2 MaxBraid™ Sutures with Needles (Pkg. 2)

## 5.5mm LactoScrew® LS Suture Anchors

<b>905490</b>	Two #2 MaxBraid™ Sutures (Pkg. 2)
<b>905496</b>	Two #2 MaxBraid™ Sutures
<b>CP905596</b>	Two #2 Osteoprene™ Sutures with Needles
<b>CP905597</b>	Two #2 Osteoprene™ Suture

## 6.8mm LactoScrew® Suture Anchors

<b>905580</b>	Two #2 MaxBraid™ Sutures
<b>905581</b>	Two #2 MaxBraid™ Sutures with Cutting Needles

## 3.0mm Ti-Screw Suture Anchors

<b>902554</b>	One #2 Blue Suture
---------------	--------------------

## 5.0mm Ti-Screw Suture Anchors

<b>902556</b>	Two #2 Blue Suture
---------------	--------------------

## 3.0mm Ti-Screw EasySlide™ Suture Anchors

<b>902574</b>	One #2 Blue Suture
<b>902578</b>	One #2 Blue Suture with Needles
<b>902482</b>	Two 2-0 MaxBraid™ Sutures
<b>902482</b>	Two 2-0 MaxBraid™ Sutures
<b>902569</b>	One #2 MaxBraid™ Sutures with Needles
<b>902570</b>	One #2 MaxBraid™ Suture

## 5.0mm Ti-Screw EasySlide™ Suture Anchors

<b>902576</b>	Two #2 Sutures
<b>902579</b>	Two #2 Sutures with Needles
<b>902571</b>	Two #2 MaxBraid™ Sutures with Needles
<b>902572</b>	Two #2 MaxBraid™ Sutures
<b>902596</b>	Two #2 MaxBraid™ Sutures with Needles

## 6.5mm Ti-Screw EasySlide™ Suture Anchors

<b>902580</b>	Two #2 Sutures with Needles
<b>902573</b>	Two #2 MaxBraid™ Sutures
<b>902587</b>	Two #2 MaxBraid™ Sutures with Needles

## 5.0mm Ti-Screw SP Suture Anchors

<b>902588</b>	Two #2 MaxBraid™ Sutures with Needles
<b>902591</b>	Two #2 MaxBraid™ Sutures
<b>902597</b>	Three #2 MaxBraid™ Sutures
<b>902598</b>	Three #2 MaxBraid™ Sutures with Needles

## 6.5mm Ti-Screw SP Suture Anchors

<b>902589</b>	Two #2 MaxBraid™ Sutures with Needles
<b>902592</b>	Two #2 MaxBraid™ Sutures with Needles
<b>902599</b>	Three #2 MaxBraid™ Sutures
<b>902600</b>	Three #2 MaxBraid™ Sutures with Needles

# Ordering Information

## Tray #1—Hand Instruments/ RC Retractor

### Charlotte™ Shoulder System Tray #1

903730

### Charlotte™ Shoulder System Retractor

903705 Body  
903704 Deep Blade, 1 Pair  
903706 Small Blade, 1 Pair  
903707 Large Blade, 1 Pair

### Grasper

902075 2.75mm  
902083 3.4mm

### 15° Angulated Up Low Profile Biter

902733 3.4mm

### Nordt™ Knot Tightener

901477

### MicroShears

902076 2.75mm

### Caspari™ Suture Punch

901070 Straight  
902070 Mini

### Oval Cannula

900419 50mm Small  
900421 70mm Small

### Oval Cannulated Obturator

900459 50mm Small  
900461 70mm Small

### ArthroPasser™ Retrieving Device

902800 2.75mm 15° Up  
902801 2.75mm Straight  
902802 2.75mm 45° Up

## Tray #2—Portal Placement/ Miscellaneous

### Charlotte™ Shoulder System Tray #2

903731

### Charlotte™ Shoulder System Shoulder Rasp

903709

### Guide Pin

903720 Pin  
903721 Handle

### Calibrated Probe

901010

### Charlotte™ Shoulder System Needle Driver

903708

### RC Buttress Template

905971

### RC Buttress Awl

905972

### Curved Rasp with Elevator

903719

### Reusable Cannula Obturator

900328 7mm  
900370 9mm

### Charlotte™ Shoulder System Tunnel Awl

903710 Small  
903712 Large

### Charlotte™ Shoulder System Tunnel Retriever

903714 Small  
903716 Large

## Tray #3—Anchor Instrumentation/ Suture Management

### Charlotte™ Shoulder System Tray #3

903732

### LactoScrew® Tap

905584

### Screw/Anchor Tap

905586 5.5mm

### LactoScrew® Drill

905585 3.5mm

### Screw/Anchor Awl

905588 5.5mm

### Drill Guide Obturator

905561

### Ratchet Handle

909733

### Open Knot Pusher

902803

### Suture Crochet Hook

903610

### SpeedPass™ Suture Retriever—Reusable

902804 70° Left Pigtail—Reusable  
902805 70° Right Pigtail—Reusable  
902806 45° Left Pigtail—Reusable  
902807 45° Right Pigtail—Reusable  
902808 70° Straight Up—Reusable

## Miscellaneous

### Charlotte™ Shoulder System Tray #5

903734

### Charlotte™ Shoulder System Tray Lid

903740

## Caspari™ Suture Punch

### Caspari™ Suture Punch

901070 Straight  
901070HP Straight—High Performance

### Caspari™ Mini Suture Punch

902070 Straight  
902070HP Straight—High Performance  
902071 Left  
902071HP Left—High Performance  
902072 Right  
902072HP Right—High Performance  
902073 15° Up  
902073HP 15° Up—High Performance

### Caspari™ Suture Punch/ Hand Instrument Case

902095

## Nordt™ Knot Tightener

### Nordt™ Knot Tightener

901477  
901477HP High Performance

## Ti-Screw Instruments

### 1.5mm Drill Bit

# Ordering Information

902562

## Ti-Screw Anchor Tap

---

902565 5.0mm

902566 3.0mm

## Additional LactoScrew® Instruments

### LactoScrew® Anchor Tap

---

905587 6.8mm

905497 2.8mm

### LactoScrew® Anchor Drill

---

905492 2.8mm

*Charlotte,™ MaxBraid,™ LactoScrew,® EasySlide,™ ArthroPasser,™ Nordt™ Caspari™ and MaxCutter™ are trademarks of Arthrotek, Inc.*

*LactoSorb® is a trademark of Biomet Manufacturing Corp.*

*Dyneema® is a trademark of DSM IP Assets BV.*

*Osteoprene™ is a trademark of Poly-Med, Inc.*

This material is intended for the Arthrotek Sales Force and surgeons only.  
It is not intended to be redistributed without the express written consent of Arthrotek.



P.O. Box 587, Warsaw, IN 46581-0587 • 800.348.9500 ext. 1501 • ©2006 Arthrotek, Inc. All Rights Reserved  
web site: [www.arthrotek.com](http://www.arthrotek.com) • eMail: [arthrotek@arthrotek.com](mailto:arthrotek@arthrotek.com)  
Form No. Y-BMT-960/033106/K

I N V E N T I N G   T H E   F U T U R E   O F   A R T H R O S C O P Y

