

SLAP LESION REPAIR

Using the 3.5mm LactoScrew™ Suture Anchor

Surgical Technique

Diagnostic Arthroscopy

- Standard posterior portal (Figure 1)
- Standard anterior rotator interval portal (using spinal needle to assure proper angle to address anterior/superior labrum)
- Superior labral probing: assess for SLAP lesion
 - **Signs of injury (labral fraying, exposed bone)**
 - **Superior displacement of labrum > 5mm with probing**
 - **Inferior displacement**
- Biceps tendon: assess for partial-thickness tearing/instability
 - **Origin**
 - **Intraarticular component**
 - **Reduce tendon in the joint**
 - **Rule out rotator interval lesion/assess medial stability**
 - **Consider biceps tenodesis or tenolysis if significant partial-thickness tearing and/or medial instability (Figure 2)**
- Other...
 - **Rotator cuff**
 - **Articular cartilage**
 - **Ligaments**
 - **Synovium**




Figure 1

Superior Labral Lesion (SLAP)





Figure 2: Biceps Partial Tear

Superior Labral Debridement/Supraglenoid Preparation

- Contour edges
- Small, curved shaver (3.5mm) to debride superior glenoid to bone
- Small burr (3.5mm) to decorticate superior glenoid bone
 - **Avoid deep trough**
 - **Avoid glenoid articular cartilage injury**



Type II SLAP Lesion After Debridement

Surgical Technique

Establish Mid-Lateral Portal

- 2mm guide pin just off mid-lateral acromial edge (Figure 3)
- Pin transveres through myotendonous junction of supraspinatus
- Sharp 7mm cannulated obturator +7mm cannula placed over guide pin

Placement of Suture Anchor(s)

- Probe in mid-lateral portal pushing labrum medially, exposing prepared superior glenoid bone
- If two suture anchors are utilized, start with the most anterior anchor
- Drill/tap to laser mark depth through anterior portal
- LactoScrew® Suture Anchor placed to laser mark depth
- Assess stability of anchor under direct visualization
- Pull a suture limb through lateral portal to avoid entanglement

Placement of Suture Through Anterior Labrum

- Place a Mini Suture Punch through anterior portal to pierce labrum (Figure 4) (for simple suture technique)
- Transfer other limbs to anterior portal utilizing arthroscopic suture retriever
- Remove and replace anterior cannula with sutures **outside** this cannula (avoids entanglement of sutures)

If Posterior Suture Anchor is Desired

- Place LactoScrew® Suture Anchor through lateral portal
- Transfer limb of posterior suture anchor to anterior portal utilizing arthroscopic suture retriever
- Place Mini Suture Punch through lateral portal to pierce labrum (Figure 5) (for simple suture technique)
- Transfer other limb of posterior anchor back through lateral portal utilizing suture retriever




Figure 3




Figure 4





Figure 5

Knot Tying²

- Tie posterior anchor sutures first with the Nordt™ Knot Tightener
- **First throw:** suture through labrum is post
- **Second throw:** two half-hitches, same direction, around post utilizing Nordt™ Knot Tightner
- **Third throw:** half-hitch, opposite direction, alternate post (locking the knot)
- **Fourth throw:** half-hitch, opposite direction, alternate post
- **Fifth throw:** half-hitch, opposite direction, alternate post
- Cut suture ends 1–2mm above knot
- Pull second suture ends from outside anterior cannula through lateral portal
- Repeat knot tying of anterior anchor sutures
- Assess stability of repair



Nordt™ Knot Tightener



(A) Status Post Repair; (B) Assessing for Superior Stability;
(C) Assessing for Inferior Displacement

1. D'Alessandro DF, Valadie AL. Superior glenoid lesions: a diagnostic and therapeutic challenge. *J South Orthop Assoc.* 4(3): 214-27, 1995.

2. Connor PM, Milia M. Arthroscopic knot tying techniques: a critical assessment. *J Arthroscopy*, Submitted for publication, 2000.

Developed in conjunction with Patrick M. Connor, M.D. and Donald F. D'Alessandro, M.D., Charlotte, North Carolina.

Charlotte,™ LactoScrew™ and Nordt™ are trademarks of Arthrotek, Inc.

This brochure is presented to demonstrate the surgical technique utilized by Patrick M. Connor, M.D. and Donald F. D'Alessandro, M.D. Arthrotek, as the manufacturer of this device, does not practice medicine and does not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs any procedure is responsible for determining and utilizing the appropriate techniques for such procedure for each individual patient. Arthrotek is not responsible for selection of the appropriate surgical technique to be utilized for an individual patient.

This material is intended for the Arthrotek Sales Force and surgeons only. It is not intended to be redistributed without the express written consent of Arthrotek.

Arthrotek, Inc.
A Wholly Owned Subsidiary of Biomet, Inc.
56 East Bell Drive
Warsaw, Indiana 46582 USA

01-50-1072
Date: 02/05

ARTHROTEK® SOFT TISSUE ANCHORING DEVICES
ATTENTION OPERATING SURGEON

DESCRIPTION

The Arthrotek® Soft Tissue Anchoring Devices are resorbable repair devices used to attach soft tissue to bone. LactoSorb® Soft Tissue Screw and Washer and MicroMax™ Suture Anchor are used with or without a suture. LactoSorb® L-15 Screw Anchors consist of a screw and head design and are used with a suture. The devices are implanted into a predrilled bone hole and are made of a resorbable copolymer, a polyester derivative of lactic acid and glycolic acid. Polylactic/polyglycolic acid copolymer degrades and resorbs in vivo by hydrolysis to lactic and glycolic acids, which are then metabolized by the body.

MATERIALS

Poly-L-Lactic Acid/Polyglycolic Acid
Ultra-High Molecular Weight Polyethylene (UHMWPE)
Polyester
Polypropylene

INDICATIONS

1. LactoScrew® L-15 Screw Anchor (85% PLLA/15% PGA):

Shoulder

Bankart repair
SLAP lesion repair
Acromio-clavicular separation
Rotator cuff repair
Capsule repair or capsulolabral reconstruction
Biceps tenodesis
Deltoid repair

Wrist/Hand

Scapholunate ligament reconstruction
Ulnar/radial collateral ligament reconstruction

Ankle/Foot

Lateral stabilization
Medial stabilization
Achilles tendon repair/reconstruction
Hallux valgus reconstruction
Mid- and forefoot reconstruction

Elbow

Tennis elbow repair
Ulnar or radial collateral ligament reconstruction
Biceps tendon reconstruction

Knee

Medial collateral ligament repair
Lateral collateral ligament repair
Posterior oblique repair
Joint capsule closure
Iliotibial band tenodesis

Patellar ligament/tendon repair

2. LactoSorb® L-15 Screw and Washer (85% PLLA/15% PGA) and MicroMax™ Suture Anchor:

Shoulder Indications

Bankart Repair
SLAP Lesion Repair
Acromio-clavicular Separation Repair
Rotator Cuff Repair
Capsule Repair and Capsulolabral Reconstruction
Biceps Tenodesis
Deltoid Repair
Wrist Indications
Scapholunate ligament reconstruction

Elbow Indications

Tennis Elbow Repair
Biceps Tendon Reattachment
Medial and Lateral Repairs
Ulnar or Radial Collateral Ligament Reconstruction

Knee Indications

Extra-Capsular Repair
Medial Collateral Ligament Repair
Lateral Collateral Ligament Repair
Posterior Oblique Ligament Repair
Joint Capsule Closure
Iliotibial Band Tenodesis Reconstruction
Patellar Ligament/Tendon Repair
Vastus Medialis Obliquus (VMO) Muscle Advancement

LactoSorb® L-15 Screw and Washer and MicroMax™ Suture Anchor are preloaded with suture for use at the discretion of the physician.

CONTRAINDICATIONS

1. Active infection.
2. Patients with mental or neurologic conditions who are unwilling or incapable of following postoperative care instructions.

3. Patient conditions including: blood supply limitations, insufficient quantity or quality of bone, or latent infections.
4. Pathologic soft tissue conditions, which would prevent secure fixation.

WARNINGS

Arthrotek® internal fixation devices provide the surgeon with a means to aid in the management of soft tissue to bone reattachment procedures. While these devices are generally successful in attaining these goals, they cannot be expected to replace normal healthy soft tissue or withstand the stress placed upon the device by full or partial weight bearing or load bearing, particularly in the presence of incomplete healing. Therefore, it is important that immobilization (use of external support, sling, etc.) of the treatment site be maintained until healing has occurred. Surgical implants are subject to repeated stresses in use, which can result in fracture or damage to the implant. Factors such as the patient's activity level and adherence to weight bearing or load bearing instructions have an effect on the service life of the implant. The surgeon must be thoroughly knowledgeable not only in the medical and surgical aspects of the implant, but also must be aware of the mechanical and polymeric aspects of the surgical implants.

1. Correct selection of the implant is extremely important. The potential for success in soft tissue to bone fixation is increased by the selection of the proper type of implant. While proper selection can help minimize risks, the device is not designed to withstand the unsupported stress of full weight bearing, load bearing or excessive activity.
2. The implants can loosen or be damaged when subjected to increased loading associated with inadequate healing. If healing is delayed, or does not occur, the implant or the procedure may fail. Loads produced by weight bearing and activity levels may dictate the longevity of the implant.
3. Inadequate fixation at the time of surgery can increase the risk of loosening and migration of the device or tissue supported by the device. Sufficient bone quantity and quality are important for adequate fixation and success of the procedure. Bone quality must be assessed at the time of surgery. Adequate fixation in diseased bone may be more difficult. Patients with poor quality bone, such as osteoporotic bone, are at greater risk of device loosening and procedure failure.
4. Care is to be taken to assure adequate soft tissue fixation at the time of surgery. Failure to achieve adequate fixation or improper positioning or placement of the device can contribute to a subsequent undesirable result.
5. The use of appropriate immobilization and postoperative management is indicated as part of the treatment until healing has occurred.
6. Correct handling of implants is extremely important. Do not modify implants. Do not notch or bend implants. Notches or scratches put in the implant during the course of surgery may contribute to breakage. Intraoperative fracture of devices can occur if excessive force (torque) is applied while seating.
7. DO NOT USE if there is loss of sterility of the device.
8. Discard and DO NOT USE opened or damaged devices, and use only devices that are packaged in unopened or undamaged containers.
9. Ensure contact of tissue to bone when implanting. DO NOT OVERTIGHTEN the screw. Structural damage to the tissue and implant may occur if the screw is overtightened.
10. Adequately instruct the patient. Postoperative care is important. The patient's ability and willingness to follow instructions is one of the most important aspects of successful soft tissue management. Patients affected with senility, mental illness, alcoholism, and drug abuse may be at a higher risk of device or procedure failure. These patients may ignore instructions and activity restrictions. The patient is to be instructed in the use of external supports that are intended to immobilize the repair site and limit weight bearing or load bearing. The patient is to be made fully aware and warned that the device does not replace normal healthy tissue, and that the device can break, bend or be damaged as a result of stress, activity, load bearing, or weight bearing. The patient is to be made aware and warned of general surgical risks, possible adverse effects, and to follow the instructions of the treating physician. The patient is to be advised of the need for regular postoperative follow-up examination as long as the device remains implanted.
11. MicroMax™ Suture Anchor—Loss of bone fixation may occur if flanged wings are not properly deployed.

PRECAUTIONS

Instruments are available to aid in the accurate implantation of internal fixation devices. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments, which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Arthrotek recommends that all instruments be regularly inspected for wear and disfigurement.

If device contains MaxBraid™ suture, refer to manufacturer package insert for further information.

POSSIBLE ADVERSE EFFECTS

1. Infection can lead to failure of the procedure.
2. Neurovascular injuries can occur due to surgical trauma.
3. Bending, fracture, loosening, rubbing, and migration of the implant may occur as a result of excessive activity, trauma, or load bearing.
4. Implantation of foreign materials can result in an inflammatory response or allergic reaction.
5. Inadequate healing.
6. Pain, discomfort, or abnormal sensation due to the presence of the device.
7. Necrosis of the bone or tissue.

STERILITY

Arthrotek® resorbable implants are sterilized by exposure to Ethylene Oxide (ETO) Gas. Do not resterilize. Do not use past expiration date.

STORE AT OR BELOW ROOM TEMPERATURE. DO NOT EXPOSE PRODUCT TO TEMPERATURES GREATER THAN 120°F OR 49°C.

Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.
Authorized Representative: Biomet U.K., Ltd.
Waterton Industrial Estates
Bridgend, South Wales
CF31 3XA U.K.

CE 0086

The information contained in this package insert was current on the date this brochure was printed. However, the package insert may have been revised after that date. To obtain a current package insert, please contact Arthrotek at the contact information provided herein.

ARTHROTEK®
A BIOMET COMPANY

P.O. Box 587, Warsaw, IN 46581-0587 • 800.348.9500 ext. 1501 • ©2002, 2006 Arthrotek, Inc. All Rights Reserved
web site: www.arthrotek.com • eMail: arthrotek@arthrotek.com

Form No. Y-BMT-963/033106/K

I N V E N T I N G T H E F U T U R E O F A R T H R O S C O P Y

