
ENCORE® Shoulder System
ENCORE® Schulter-System
Système de humérale ENCORE®
Sistema de humeral ENCORE®
Sistema de protesi della spalla ENCORE®

CE 0123

MDSS, Burckhardtstrasse 1
D-30163 Hannover, Germany
Telefon 49 511 6262 8630

0400-0061 F

E

1. Product Handling

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage which may compromise sterility. If packaging has been opened or damaged, contact manufacturer's representative. When unpacking the implant, verify the labeling for correct Ref. No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects which may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original FOUNDATION® product belonging to the same system.

2. Product Description and Implant Materials

The FOUNDATION® Shoulder System is a semi-constrained total shoulder replacement consisting of humeral and glenoid components for cemented or cementless application. The humeral stem prosthesis is manufactured from titanium alloy (wrought Ti6Al4V per ASTM F136 or cast titanium per ASTM1108) and is available with plasma spray, roughened surface or porous coating on the proximal body. Humeral head components are available in CoCr (ASTM F799). All polyethylene glenoid components are manufactured from UHMWPE (ASTM F648, ISO Standard 5834/1+2).

3. Indications

Total joint replacement is indicated for patients suffering from disability due to:

- noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis of the natural humeral head and/or glenoid;
- rheumatoid arthritis;
- correction of functional deformity;
- humeral fracture.

This device may also be indicated in the salvage of previously failed surgical attempts.

While total shoulder replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

4. Contraindications

Total joint replacement is contraindicated where there is:

- infection or sepsis;
- insufficient bone quality which may affect the stability of the implant;
- muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- alcoholism or other addictions;
- materials sensitivity;
- loss of ligamentous structures;
- high levels of physical activity (e.g. competitive sports, heavy physical labor)

5. Precautions and Warnings

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Only FOUNDATION® Shoulder System instruments and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

6. Preoperative Planning and Postoperative Care

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the FOUNDATION® Shoulder system are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his/her activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

7. Adverse Effects

1. Accelerated wear of the polyethylene articulating surfaces have been reported following total shoulder replacement. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Accelerated wear shortens the useful life of the prosthesis, and leads to early revision surgery to replace the worn prosthetic components.
2. Metallosis and osteolysis may be implicated from wear debris associated with the use of orthopedic implants.
3. Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage occurs more frequently, possibly the result of surgical

trauma.

4. Metal sensitivity reactions in patients following joint replacement have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. In some cases, wear debris can initiate the process of histiocytic granuloma formation and consequent osteolysis and loosening of the implant.
5. Dislocation and subluxation of implant components can result from improper positioning of the components. Muscle and fibrous tissue laxity can also contribute to these conditions.
6. Implants can loosen or migrate due to trauma or loss of fixation.
7. Infection can lead to failure of the joint replacement.
8. While rare, fatigue fracture of the implant can occur as a result of strenuous activity, improper alignment, or duration of service.
9. Fracture of the humerus can occur while press-fitting (seating) the humeral stem into the prepared humeral canal.
10. Allergic reactions.

Intraoperative and early postoperative complications can include:

- 1) humeral perforation, or fracture;
- 2) humeral fracture can occur while seating the device;
- 3) damage to blood vessels;
- 4) temporary or permanent nerve damage resulting in pain or numbness of the affected limb;
- 5) undesirable shortening or lengthening of the limb;
- 6) traumatic arthrosis of the shoulder from intraoperative positioning of the extremity;
- 7) cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- 8) hematoma;
- 9) delayed wound healing; and,
- 10) infection.

Late postoperative complications can include:

- 1) avulsion as a result of excess muscular weakening;
- 2) non-union due to inadequate reattachment and/or early weight bearing;
- 3) aggravated problems of other joints of the affected limb or muscle deficiencies;
- 4) humeral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- 5) periarticular calcification or ossification, with or without impediment to joint mobility;
- 6) inadequate range of motion due to improper selection or positioning of components, by impingement, and calcification.

8. Sterilization

Unless opened or damaged, Encore® implants are supplied sterile in multiple pouches or barrier blister trays. Check all packaging for punctures or other damage. If packaging is opened or damaged, contact manufacturer or manufacturer's representative for instructions.

Sterilization of implants is by gamma radiation at the minimum dose of 25 kGy to achieve a Sterility Assurance Level (SAL) of 10^{-6} . Implants are single-use devices. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. However, should the original sterile package be inadvertently opened or compromised before implantation, the metallic devices may be resterilized prior to use following the guidelines listed below.

Do not resterilize an implant or component that has been in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from porous coating and storage of the implant or component should be avoided. The manufacturer and distributor assume no responsibility for the cleaning and/or resterilization of implants, components, or reusable instruments performed by the individual or hospital.

Encore instruments and instrument cases are generally composed of titanium, stainless steel, aluminum, and/or polymeric materials. The cases may be multi-layered with various inserts to hold surgical instrumentation in place during handling and storage. The inserts may consist of trays, holders, and silicone mats. The instrument cases will allow sterilization of the contents to occur in a steam autoclave utilizing the cleaning, sterilization, and drying cycle that has been validated and listed below. Instrument cases do not provide a sterile barrier and must be used in conjunction with sterilization wrap to maintain sterility. Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines listed below.

User Resterilization Guidelines for Implants:

Porous coated or non-coated metallic implants or components, which have not been previously used or implanted, can be resterilized using the following steam sterilization cycle. All devices should be disassembled prior to resterilization and care should be taken to protect implant or component from mechanical damage.

WARNING: DO NOT resterilize UHMWPE (ultra-high molecular weight polyethylene) implants, PMMA (polymethylmethacrylate) spacers, HA (hydroxylapatite) coated implants, and ceramic implants.

<u>Steam Resterilization</u> for metallic implants and components ONLY!	Remove from supplied packaging and wrap in protective sterilization wrap according to AAMI / AORN guidelines or place into appropriate case configuration. Pre-Vacuum Autoclave (HI-VAC): 270-272° F (132-134° C), 6-minute exposure time, with 4 pulses and a 5-minute dry time shall achieve a SAL of 10^{-6} . Gravity Displacement Autoclave: 270-272° F (132-134° C), 30-minute exposure time, with a 30-minute dry time shall achieve a SAL of 10^{-6} .
---	--

Encore has validated the above steam sterilization cycles and has data on file. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques. Proper validation of the autoclave is essential to ensure proper sterilization temperatures and cycle times. **NOTE: Encore does not recommend Flash or Chemical Sterilization.**

User Cleaning and Resterilization Guidelines for Reusable Surgical Instruments:

Keep instruments moist and do not allow blood and/or bodily fluids to dry on the instruments. The decontamination process should begin immediately after completion of the surgical procedure.

All instruments should be positioned to allow sterilant to come in contact with all surfaces. All jointed instruments should be in the open or unlocked position with ratchets not engaged. Instruments composed of more than one part, with sliding pieces, screws, or removable parts should be disassembled, if possible.

1. **Decontamination:** Saturate the surface completely with full strength disinfectant/cleaner* (e.g. Cavicide) and allow to remain in contact with devices for 5 minutes.
2. **Pre-Cleaning:** Remove gross contaminants by immersing the devices in room temperature neutral pH enzymatic cleaner* (e.g. Metrizyme) and disassemble instruments, if suitable. Scrub with the appropriate soft bristle brush until visibly clean.
3. **Washing:** Immerse devices in the ultrasonic washer/cleaner with room temperature neutral pH enzymatic cleaner* (e.g. Metrizyme) and sonicate for 10 minutes. For ultrasonic cleaning follow the manufacturer's specifications for suggested water level and concentration. When using mechanical washers, make sure the instruments are secured in place, and do not touch or overlap.
4. **Rinsing:** Thoroughly rinse the devices with deionized or distilled water. For example, a minimum of 2 minutes, three (3) times.
5. **Drying:** Allow devices to air dry a minimum of 20 minutes prior to inspection and sterilization preparation. Instruments must be thoroughly dried to remove residual moisture before they are stored.
6. **Preparation and Assembly:** After cleaning/disinfection, the disassembled instruments should be reassembled and visually inspected. Check for misalignment, burrs, bent, or fractured tips. Mechanically test the working parts to verify that each instrument functions correctly. Place instruments into appropriate configuration within instrument case and wrap with protective sterilization wrap according to AAMI / AORN guidelines.

* Do not use high acidic (pH <4) or high alkaline (pH >10) products for disinfection or cleaning, since these can corrode metal, cause discoloration or stress fractures. Encore has qualified the above cleaning method with the provided solution examples, for a 3 Spore Log Reduction (SLR). Other cleaning/disinfection methods may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques.

7. **Sterilization with a Pre-Vacuum Autoclave (HI-VAC):**
270-272° F (132-134° C), 16-minute exposure time, with 4 pulses and a 30-minute dry time.
Sterilization with a Gravity Displacement Autoclave:
270-272° F (132-134° C), 30-minute exposure time, with a 30-minute dry time.

Encore has validated the above steam sterilization cycles achieve a SAL of 10⁻⁶ with components loosened and has data on file. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques. Proper validation of the autoclave is essential to ensure proper sterilization temperatures and cycle times. **NOTE: Encore does not recommend Flash Sterilization within instrument cases or Chemical Steriliza**

Warning: Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating. For further information regarding the use of Encore Shoulder Systems, contact your distributor. Encore Medical, L.P. (Made in USA)

9. Trademarks and patents

FOUNDATION® is a trademark of Encore Medical, L.P.

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