

Trace minerals for Methanococcus medium. Nitrogen-free. 1,000X.

	<u>g/100ml</u>
Na <sub>3</sub> Citrate·2H <sub>2</sub> O	2.1
Adjust pH to 6.5	
MnSO <sub>4</sub> ·2H <sub>2</sub> O	0.5
CoSO <sub>4</sub> or CoCl <sub>2</sub> (·6H <sub>2</sub> O)	0.1
ZnSO <sub>4</sub> (·7H <sub>2</sub> O)	0.1
CuSO <sub>4</sub> ·5H <sub>2</sub> O	0.01
AlK(SO <sub>4</sub> ) <sub>2</sub>	0.01
H <sub>3</sub> BO <sub>4</sub>	0.01
Na <sub>2</sub> MoO <sub>4</sub> ·2H <sub>2</sub> O	0.1
NiCl <sub>2</sub> ·6H <sub>2</sub> O	0.025
Na <sub>2</sub> SeO <sub>3</sub>	0.2
V(III)Cl	0.01
Na <sub>2</sub> WO <sub>4</sub> ·2H <sub>2</sub> O	0.0033