

Negative staining protocol

1. Staining solution preparation (refer to the protocols below)
2. Glow-discharge EM grids coated with continuous carbon film (the grids should be used within 30 min after the glow-discharging). Note: usually 20mA* 20s will be used.
3. Apply 6 μ l of protein sample solution to the EM grid and wait for 1 minute.
4. Blot the grid from the edge with a piece of filter paper (briefly touch the rim of the grid with a piece of filter paper, allow the solution to flow to the filter paper). Caution: don't let the grid to be completely dry. It is better to have a thin layer of solution.
5. Apply 3-6 μ l staining solution and wait for 1-2 minutes.
6. Blot with a piece of filter paper (briefly touch the rim of the grid with a filter paper, allow the solution to go into the filter paper completely).
7. Completely dry the grid in the hood for 10 minutes.
8. Store the stained grid in a grid box and label properly.

Staining solution preparation

1. 2% Uranyl Acetate

- 1.1. Weight 0.2 g of uranyl acetate
- 1.2. Stir to dissolve the uranyl acetate particles in 10 ml ddH₂O (can take hours) in the hood.
- 1.3. Filter solution with a 0.22 μ m filter using a syringe
- 1.4. Aliquot into eppendoff tubes and wrap with aluminum foil
- 1.5. Store at room temperature in the dark
- 1.6. Filter again before use after a few days

2. 2% Uranyl formate

- 2.1. Boil 10ml ddH₂O and let it cool to room temperature
- 2.2. Weight 20mg uranyl formate in Eppendorf tube.
- 2.3. Add 1ml water to the tube.
- 2.4. Vortex vigorously to dissolve the yellow powder
- 2.5. Add 10ul 1N NaOH to the cap of the tube.
- 2.6. Close the cap and vortex vigorously
- 2.7. Filter this solution through a 0.02 μ m filter using a syringe

Here is another protocol used at LBNL to make 1% UF

- 2.1) Put 1 g Uranyl Formate (UF) powder into a glass bottle containing 100 ml deionized water.
- 2.2) Stir solution overnight at room temperature under a dark condition, in which the water bottle was wrapped with aluminum foil to prevent light.

- 2.3) Aluminum foil wrap a 5 ml syringe using a 0.2 μm (pore size) filter, and filter the solution gently through the syringe into an aluminum foil covered tube.
- 2.4) Aluminum-foil wrap a 1 ml syringe using a 0.02 μm (pore size) filter, then filter the solution followed by aliquot the filtered solution into 2 ml vials covered by aluminum foil protecting the stain from light.
- 2.5) Immediately after aliquot, freeze the vials by plunging the vials into liquid nitrogen container using long handled forceps.
- 2.6) Transfer the frozen vials into a $-80\text{ }^{\circ}\text{C}$ freezer for storage for future usage.
- 2.7) Thaw a vial of 1% UF solution in a water bath set at room temperature. Make sure that the aluminum foil remains wrapped around the vial to prevent exposure to light.
- 2.8) Once the UF is completely thawed, using an aluminum foil wrapped 1 ml syringe mounted with a 0.02 μm filter to filter the UF solution right before usage. Save the filtered solution into a new aluminum foil wrapped vial, and then place the vial into a cap covered ice box to maintain a dark environment.

3. **2% Sodium Silicotungstate (SST)**

- 3.1. Weight 0.2 g sodium silicotungstate and put into 10 ml ddH₂O
- 3.2. Stir to dissolve sodium silicotungstate particles (it can take hours to overnight)
- 3.3. Adjust pH using couple drops of NaOH. (Adjust PH to 7.0?)
- 3.4. Filter solution with a 0.22 μm filter using a syringe
- 3.5. Aliquot into eppendoff tubes and wrap with aluminum foil
- 3.6. Store at room temperature in the dark

4. **5% Ammonium molybdate (AM)**

- 4.1. Weight 0.05 g ammonium molybdate and place in a 1.5 ml Eppendorf tube
- 4.2. Add 0.8 ml distilled water and dissolve all the ammonium molybdate
- 4.3. Add a few drops of 10 N NaOH and adjust the solution pH to 7.0 (using pH paper)
- 4.4. Bring the final volume to 1.0 ml
- 4.5. Filter this solution through a 0.22 μm filter using a syringe.