



4. Using scissors, remove half of a young leaf from the first plant to be genotyped (typically, plant A for any given column). Using a **clean** 0.5mm Harris Uni-Core punch and a **clean portion** of cutting mat, collect one leaf disk and eject it into Row A under the appropriate column. Collect a second leaf disk and eject it into Row B. Verify that the disks are in the solution and not stuck to the side of the well, or missing.

5. Place the used punch in a rack in a container with 2% bleach of sufficient depth to just cover the steel tips of the punch.

6. Use a clean punch and a clean portion of cutting mat for each additional sample.

7. When all leaf disks have been collected, put the 384-well plate into a thermocycler. Cover the plate with a FRESH single sheet of plate sealer (clear vinyl from a fabric store, cut to size, is the cheapest). Close the lid on the thermocycler.

8. Run the “phire” PCR program:

1 cycle of 98C x 5min

40 cycles of 98C x 5sec, 58C x 5sec, 72C x 20sec

1 cycle of 72C x 1min

The whole program takes just over an hour to complete. After PCR the plate may be stored at 4C for several days if necessary.